

Cell Meter 细菌活性检测试剂盒

Table 1 Contents and storage

Material	Amount	Storage	Stability
MycoLig™ Green	1 vial (200 µL)	<ul style="list-style-type: none"> • ≤ -20°C • Desiccate • Protect from light 	
Propidium iodide	1 vial (200 µL)		
Spectral characteristic of the fluorescent probe: Ex~510, Em~530			

Introduction

Mycolight™ Bacterial Viability Assay Kit provides two-color fluorescence assay of bacterial viability in both gram-positive and negative bacterial cell. The kit utilizes the mixture of our green fluorescent nucleic acid stain MycoLight™ Green and the red-fluorescent nucleic acid stain propidium iodide. When used alone, the MycoLight™ Green stain generally labels all bacteria (live and dead) in a population. In contrast, propidium iodide penetrates only bacteria with damaged membranes, causing a reduction in the MycoLight™ Green stain fluorescence when both dyes are present. Thus, with an appropriate mixture of the MycoLight™ Green and propidium iodide stains, live bacteria with intact cell membranes emits green fluorescence, whereas dead or dying bacteria with damaged membranes gives red fluorescence. The Mycolight™ Bacterial Viability Assay Kit is a robust tool for monitoring the viability of bacterial populations as a function of the membrane integrity of the cell. Stained cells can be monitored fluorimetrically at 510-530 nm (FITC filter) and 600-660 nm (Texas red filter) with excitation at 488 nm, the most common excitation light source.

Guidelines for Use

Important Thaw kit components at room temperature and centrifuge briefly before starting

your experiment.

The Kit has been tested at logarithmically growing cultures of the following bacterial species: *Bacillus cereus*, *B. subtilis*, *Clostridium perfringens*, *Escherichia coli*, *Klebsiella pneumoniae*, *Micrococcus luteus*, *Mycobacterium phlei*, *Pseudomonas aeruginosa*, *P. syringae*, *Salmonella oranienburg*, *Serratia marcescens*, *Shigella sonnei*, *Staphylococcus aureus* and *Streptococcus pyogenes*. *Agrobacterium tumefaciens*, *Edwardsiella ictaluri*, *Eurioplasma eurilytica*, *Lactobacillus sp.*, *Mycoplasma hominus*, *Propionibacterium sp.*, *Proteus mirabilis* and *Zymomonas sp.*

The following is the recommended protocol for bacterial staining. The protocol only provides a guideline and should be modified according to the specific needs.

experiment procedure

1. Grow bacteria in any appropriate medium. Best results for healthy bacteria are obtained from log-phase cultures. Dilute the bacterial culture to ~ 10⁶ to 10⁸ cells per mL in 0.85% NaCl or appropriate buffer. Prepare sufficient suspension to provide 500 µL per test for flow cytometry or 100 µL per test for 96-well plate.

Remove traces of growth medium before staining bacteria. A single wash step is usually sufficient to remove significant traces of interfering

BIOLUMINOR

media components from the bacterial suspension. Phosphate wash buffers are not recommended because they appear to decrease staining efficiency.

2. Add 4 μL of the dye working solution (250X) to each mL of the bacterial suspension. Mix well and incubate at room temperature for 15 minutes. Protect from light.

3. The stained bacterial cells can be analyzed by a fluorescence microscope, fluorescent microplate reader or flow cytometry.

4. The fluorescence from both live and dead bacteria may be viewed simultaneously with any standard fluorescein long pass filter set. Alternatively, the live (green fluorescent) and dead (red fluorescent) cells may be viewed separately with FITC and Texas Red filter sets.

Fluorescence Data

Instrument: Fluorescence microscope

Excitation: 510/600 nm

Emission: 530/660 nm

Recommended plate: Black wall/clear bottom

Instrument specification(s): FITC/Texas Red filter sets

Instrument: Flow cytometer

Excitation: 488 nm laser

Emission: 530/30 nm, 610/20 nm filter

Instrument specification(s): FITC, PE-Texas Red channel